

Extras of introduction to Single-cell RNA-seq

Konstantin Zaitsev, ITMO University

Masters Program, Feb 12th, 2020

Latest things: single-nuclei RNA-seq


- ✓ Not much different in sequencing and preparation
- ✓ Can be done with frozen samples
- ✓ People do that a lot with brain biobanks

Latest things: 5' scRNA-seq allows TCR/IG capture


- ✓ 5' scRNA-seq is different from 3' scRNA-seq
- ✓ 5' scRNA-seq allows us to get TCR and IG repertoires of the cells

3' scRNA-seq

Inside individual GEMs




Pooled cDNA processed in bulk




5' scRNA-seq

Inside individual GEMs



Pooled cDNA processed in bulk



Latest things: scRNA-seq + surface protein expression

BRIEF COMMUNICATIONS

nature
biotechnology

rights reserved.

Multiplexed quantification of proteins and transcripts in single cells

Vanessa M Peterson^{1,5}, Kelvin Xi Zhang^{2,5},
Namit Kumar¹, Jerelyn Wong³, Lixia Li¹, Douglas
C Wilson³, Renee Moore⁴, Terrill K McClanahan³,

the standard 10x Genomics single-cell (sc)RNA-seq platform³, which is a droplet-based system designed for 3' digital counting of mRNA in thousands of single cells.


REAP-seq leverages the DNA polymerase activity of the reverse transcriptase to simultaneously extend the primed AbB with the poly(dT) cell barcode and synthesize complementary DNA from mRNA in the same reaction. Exonuclease I is then used to degrade any excess unbound single-stranded oligonucleotides from the protein double-stranded (ds) DNA (~155 bp) products to prevent crosstalk between AbBs and cell barcodes from different cells (Supplementary Fig. 4). Dextran sulfate was added to AbB labeling buffer to reduce non-specific binding of nega-

Brief Communication

Simultaneous epitope and transcriptome measurement in single cells


Marlon Stoeckius[✉], Christoph Hafemeister, William Stephenson, Brian Houck-Loomis, Pratip K Chattopadhyay, Harold Swerdlow, Rahul Satija & Peter Smibert

Latest things: scRNA-seq + surface protein expression




Latest things: scRNA-seq + surface protein expression

Differentiation of PTPRC Expressing PBMCs via Detection of CD45RA and CD45RO Isoforms



Latest things: scRNA-seq + surface protein expression


Similar Cell Populations are Revealed Comparing Flow Cytometry and Feature Barcoding Technology Data




Latest things: scRNA-seq + pooled CRISPR screening

High-throughput Perturbation Studies Enabled via Feature Barcoding Technology

A.




B.




A. K562 cells expressing CRISPR/Cas9 were transduced with a pool of 18 gRNA's (8 targeting: 4 genes, 2 targets per gene and 10 non-targeting guides) before 10260 cells were profiled with the Chromium Single Cell Gene Expression Solution. **B.** 60% of the cells profiled were found to contain a single gRNA, 30% of cells were found to contain 2 - 6 gRNAs and the remaining 10% of cells did not contain a gRNA.


Latest things: all the things combined



Latest things: all the things combined





Latest things: all the things combined



Latest things: spatial transcriptomics

<https://science.sciencemag.org/content/353/6294/78>



Latest things: spatial transcriptomics

<https://science.sciencemag.org/content/353/6294/78>

